

pET-21a-d(+) Vectors



	Cat. No.
pET-21a DNA	69740-3
pET-21b DNA	69741-3
pET-21c DNA	69742-3
pET-21d DNA	69743-3

pET-21a(+) sequence landmarks

T7 promoter	311-327
T7 transcription start	310
T7•Tag coding sequence	207-239
Multiple cloning sites (<i>Bam</i> H I - <i>Xho</i> I)	158-203
His•Tag coding sequence	140-157
T7 terminator	26-72
<i>lac</i> I coding sequence	714-1793
pBR322 origin	3227
<i>bla</i> coding sequence	3988-4845
f1 origin	4977-5432

The maps for pET-21b(+), pET-21c(+) and pET-21d(+) are the same as pET-21a(+) (shown) with the following exceptions: pET-21b(+) is a 5442bp plasmid; subtract 1bp from each site beyond *Bam*H I at 198. pET-21c(+) is a 5441bp plasmid; subtract 2bp from each site beyond *Bam*H I at 198. pET-21d(+) is a 5440bp plasmid; the *Bam*H I site is in the same reading frame as in pET-21c(+). An *Nco* I site is substituted for the *Nde* I site with a net 1bp deletion at position 238 of pET-21c(+). As a result, *Nco* I cuts pET21d(+) at 234, and *Nhe* I cuts at 229. For the rest of the sites, subtract 3bp from each site beyond position 239 in pET-21a(+). *Nde* I does not cut pET-21d(+). Note also that *Sty* I is not unique in pET-21d(+).

The pET-21a-d(+) vectors carry an N-terminal T7•Tag[®] sequence plus an optional C-terminal His•Tag[®] sequence. These vectors differ from pET-21a-d(+) only by their selectable marker (ampicillin vs. kanamycin resistance). Unique sites are shown on the circle map. Note that the sequence is numbered by the pBR322 convention, so the T7 expression region is reversed on the circular map. The cloning/expression region of the coding strand transcribed by T7 RNA polymerase is shown below. The f1 origin is oriented so that infection with helper phage will produce virions containing single-stranded DNA that corresponds to the coding strand. Therefore, single-stranded sequencing should be performed using the T7 terminator primer (Cat. No. 69337-3).



pET-21a-d(+) cloning/expression region

pET-21a(+) Restriction Sites

Enzyme	# Sites	Locations
AccI	2	180 2935
AccII	8	831 1559 1890 2674 2815 3117 4357 5041
Acil	81	
AflIII	2	1064 3165
AluI	25	
AlwI	16	
Alw21I	9	159 190 564 1048 2159 2983 3483 4644 4729
Alw44I	4	1044 2979 3479 4725
AlwNI	1	3581
ApaI	1	1275
ApaBI	1	748
ApoI	4	192 1339 5003 5014
AvaI	1	158
Avall	7	1616 1992 2080 2171 2450 4196 4418
BamHI	1	198
BanI	9	386 407 521 984 1703 1833 1959 4006 5238
BanII	5	190 448 462 1275 5276
BbsI	4	1210 1549 1923 2283
BbvI	26	
BccI	21	
Bce83I	7	21 1878 2048 3256 3554 3795 4663
BceII	5	583 924 1551 3667 5227
BcgI	11	
BclI	1	1078
Bfal	8	70 232 277 2179 3660 3913 4248 5352
BglI	2	2128 4178
BglII	1	342
BmgI	1	1273
BpmI	5	902 1391 2025 2692 4128
Bpu10I	1	2271
Bpu1102I	1	80
BsaI	1	4119
BsaAI	2	2917 5201
BsaBI	3	341 347 2362
BsaHI	6	387 408 522 1021 1704 4595
BsaJI	6	57 501 507 1699 2137 3325
BsaWI	7	2 1383 1886 2354 3371 3518 4349
BsaXI	2	1723 5149
Bsbl	2	2881 5108
BscGI	13	
BsGI	3	915 1115 2325
Bsil	2	3338 4722
BsiEI	6	169 1849 3081 3505 4428 4577
BsII	20	
BsmAI	7	761 1166 1292 1679 2806 4119 4895
BsmBI	2	1679 2806
BsmFI	4	525 2066 2436 5416
BsoFI	46	
Bsp24I	12	
Bsp1286I	13	
BspEI	2	2 2354
BspGI	1	2691
BspLU11I	1	3165
BsrI	26	
BsrBI	4	297 3098 4899 5345
BsrDI	4	1111 1477 4119 4293
BsrFI	7	374 383 750 1962 2122 4138 5302
BssHII	1	1475
Bst1107I	1	2936
BstEII	1	1245

Enzyme	# Sites	Locations
BstXI	3	866 995 1118
BstYI	12	
Cac8I	39	
CjeI	26	
CjePI	20	
CviJI	86	
CviRI	24	
DdeI	11	
DpnI	27	
DraI	3	3924 3943 4635
DrallI	1	5201
DrdI	3	2858 3273 5156
DrdII	2	787 5206
DsaI	2	501 2137
EaeI	5	166 372 504 1738 4446
EagI	1	166
Eam1105I	1	4058
EarI	3	682 3049 4853
Ecil	4	841 3239 3385 4213
Eco47III	3	469 1970 2419
Eco57I	2	3713 4725
EcoNI	1	599
EcoO109I	3	53 497 2171
EcoRI	1	192
EcoRII	7	787 1102 1642 1699 3191 3312 3325
EcoRV	1	1514
FauI	17	
FokI	10	1110 1119 2384 2446 2524 2710 2851 4024 4205 4492
FspI	2	2146 4280
GdIII	5	166 372 504 1738 4446
HaeI	5	792 2113 3180 3191 3643
HaeII	14	
HaeIII	24	
HgaI	12	
HgiEII	2	662 3751
HhaI	45	
Hin4I	3	963 4057 4131
HincII	2	181 1570
HindIII	1	173
Hinfl	14	
HpaI	1	1570
HphI	16	
MaeII	15	
MaeIII	18	
MbolI	14	
MluI	1	1064
MmeI	3	3380 3564 5178
MnlI	26	
MseI	28	
MslI	9	1116 1404 1434 2152 2347 2738 4310 4469 4828
MspI	31	
MspA1I	9	84 1094 1664 1757 2756 2875 3507 3752 4693
MwoI	38	
NarI	4	387 408 522 1704
NciI	12	
NdeI	1	238
NgoAIV	4	374 1962 2122 5302
NheI	1	231
NlaIII	25	
NlaIV	25	
NotI	1	166
NspI	4	539 2510 2802 3169
Pfl1108I	2	1951 4076
PfIMI	1	646
PleI	9	325 613 700 1496 3059 3544 4047 5136 5144
PshAI	1	1909
Psp5II	1	2171

Enzyme	# Sites	Locations
Psp1406I	6	726 2094 2490 4284 4657 4986
PstI	1	4303
PvuI	1	4428
PvuII	3	1664 1757 2756
RcaI	3	462 3885 4893
RsaI	3	1211 2971 4538
SacI	1	190
Sall	1	179
SapI	1	3049
Sau96I	18	
Sau3AI	27	
Scal	1	4538
ScrFI	19	
SfaNI	20	
Sfcl	5	310 3430 3621 4299 5420
SgrAI	1	383
SphI	1	539
SspI	2	4862 4993
StyI	1	57
TaqI	13	
TaqII	9	972 1190 1863 3067 4406 4591 4744 4761 5105
TfiI	5	1743 2045 2215 2719 3140
ThaI	35	
TseI	26	
Tsp45I	8	1245 2073 2604 2817 2912 4314 4525 5374
Tsp509I	16	
Tth111I	1	2910
Tth111II	6	903 1596 2626 3755 3762 3794
UbaII	20	
VspI	4	325 1749 1808 4230
XbaI	1	276
XcmI	3	920 1436 1454
XhoI	1	158
XmnI	2	2723 4657

Enzymes that do not cut pET-21a(+):

AatII	AflII	AgeI	AscI	AvrII
BaeI	BseRI	BsmI	BspMI	BsrGI
Bsu36I	Clal	FseI	KpnI	MscI
MunI	NcoI	NruI	NsiI	NspV
Pacl	PmeI	PmlI	RleAI	RsrII
SacII	SexAI	SfiI	SglI	SmaI
SnaBI	SpeI	SrfI	Sse8387I	StuI
SunI	Swal			